

Primer Design Strategy

Design Methods

The two plasmids being used for developing this specific cyanobacterium are pAM4788 as the backbone vector and pCS-PET_VB190412 for the superPETase gene. For convenience's sake, they shall be referred to in this document as pAM and PET, respectively. In order to design these primers, two steps must be taken. The first step is to identify the regions of each plasmid necessary for PCR and Gibson cloning. The second step is to isolate the specific fragments of each plasmid that will be used in the PCR and cloning protocols.

Region Identification

In the PET plasmid, the region to be isolated is that of the designed superPETase gene (shown in yellow in Figure 1).




Figure 1. Graphic map of the pCS-PET_VB190412 plasmid identifying featured regions.

In the pAM plasmid, the region to be removed is the GFP gene (shown in green in Figure 2), leaving the rest of the plasmid as the backbone vector. The superPETase gene (yellow) is inserted where the GFP gene (green) was removed.




Figure 2. Graphic map of the pAM4788 plasmid identifying featured regions.

Fragment Isolation

From the entire plasmid sequence, the fragments of interest are those at the beginning and end of the identified regions. The entire sequences of both the PET plasmid and the pAM plasmid are available in the appendix. The first sequence is that of the PET plasmid with the superPETase gene text in maroon and highlighted in yellow. The second sequence is that of the pAM plasmid with the GFP gene text in green and highlighted in yellow.

PCR primer isolation

The basic requirements for a PCR primer sequence is that $n \geq 18$ and that $T_m \geq 55^\circ\text{C}$, where n is the number of base pairs and T_m is the melting temperature. A forward and reverse primer must be made for each fragment that is desired to be replicated. For the superPETase gene, the primers would be composed of the sequences of base pairs at the beginning and ending of the gene. For the pAM plasmid backbone, the primers would be composed of the sequences of base pairs just before and just after the GFP gene.

superPETase primers

Forward primer

atggcaaacacgacctattcc $n = 22$ $T_m = 55.7^\circ\text{C}$
 catgaatcacaaagtcatggcaaacacgacctattccaagccagtcgccgcgtttttctggcacagctgggggggctca < 320
 250 260 270 280 290 300 310

PET_f: ATG GCA AAC AAC GAC CTA TTC C

Reverse primer

{aactcatctcagaagaggatctgtaa} $n = 26$ $T_m = 55^\circ\text{C}$
 cgccaattgcagtgaacaaaaactcatctcagaagaggatctgtaactagataggtaatctctgctttaaagcacagaat < 1310
 1230 1240 1250 1260 1270 1280 1290 1300

PET_r: TTA CAG ATC CTC TTC TGA GAT GAG TT

pAM backbone primers

Original forward primer

aaggggtgggCGCGCCGAC n = 18 T_m = 67.1°C
 aattgtacaaataa aaggggtgggCGCGCCGACCCAGctttcttgtacaaactcggccctcaggagaaggccatcctgacg < 2560
 2490 2500 2510 2520 2530 2540 2550

pAM_f: AAG GGT GGG CGC GCC GAC

However, after consulting with NM, the primer-dimer potential became apparent as shown below:

```

5' AAGGGTGGGCGCGCCGACCA
   |||| :: :
3' ACCCAGCCGCGCGGTGGGAA
  
```

To resolve this issue, an additional 4 base pairs were added at the 3' end of the original primer

pAM_f: AAG GGT GGG CGC GCC GAC CCA G

pAM_f2

ok so we are finding some mispriming sites with pAM_f. all reactions with this primer give two bands, second product ~ 500 bp smaller.

We could possibly fix this by tweaking annealing temperature and other program features. But it is probably easiest/quickest to just fix primer.

One issue is that AM_f is in the MCS, so it is challenging to avoid primer dimer due to the large number of palindromes. Extending AM_f by 1 or 2 bp causes severe primer dimer issue (1 bp: GC bad; 2 bp AGCT terrible!). So let's extend by 3 or 4... pick 4. This is pAM_f2

```

AAGGGTGGGCGCGCCGACCCAGCTTT AM_f2
AAGGGTGGGCGCGCCGACCCAG AM_f
gggCGCGCGcgaaccagctttcttgtacaaactcggccctcaggagaaggccatcctgacggatggccttttgaagctttctagaatcgggggcccccg < 2600
2510 2520 2530 2540 2550 2560 2570 2580 2590

ggggaattctagaccgagctgggcaactgacgcctcagtgagcgggtctctcggaagcgctcgagcatggccacatgtgcttgggtgctgcgatggccagcg < 2700
2610 2620 2630 2640 2650 2660 2670 2680 2690

ttccactatcgaagctgtagggcgcgccagcgggttagccgcaaacatcgcaagccttcgttctcgatcgcgccgatcgagttgtccaaagctgatcgag < 2800
2710 2720 2730 2740 2750 2760 2770 2780 2790

ttcggctgtgtcatgctgcgggtctcgttagatgtagtgtgacctgacctgacgCGCGgtgtcacgaaatgcgcagtgattggttttgatcttctattag < 2900
2810 2820 2830 2840 2850 2860 2870 2880 2890

AAGGGTGGGCGCGCCGACCCAGCTTT AM_f2
AAGGG---TGGCGCGCGACCCAG AM_f mispriming site?? AAGGGTGGGCGCGCCGACCCAG
ggatttcgtgacagctccccagctccctgactagatcggcggttttcaagccaatgagcttcgatcgccccaagcgccagttcagcaatcgccctcattcgg < 3000
2910 2920 2930 2940 2950 2960 2970 2980 2990
  
```

pAM_f2: AAGGGTGGGCGCGCCGACCCAGCTTT

Original reverse primer

{ctttaagaaggagcccttcacc} n = 22 T_m = 55.4°C
 ccttgtttaactttaagaaggagcccttcaccatgtctaaagggtgaagaattattcactggtgtgttcccaattttggttgaattagatggtgat < 1840
 1750 1760 1770 1780 1790 1800 1810 1820 1830

pAM_r: GGT GAA GGG CTC CTT CTT AAA G

The primer-dimer potential for this sequence was very high, as shown below:

```
5'-GGTGAAGGGCTCCTTCTTAAAG-3'
      ***
3'-GAAATTCTTCCTCGGGAAGTGG-R'
```

This was resolved by adding an A at the 3' end of the original primer before converting into the complement
pAM_r: GGT GAA GGG CTC CTT CTT AAA GT

pAM Internal Primers

Running PCR amplification for the pAM backbone consist of >7k base pairs. This is too large to be efficient. To compensate for this, we divided the pAM backbone into 2 separate fragments and created primers. The primers are as follows:

Note- naming convention: 35K refers to 3.5 kbp past pAM_f

Forward primer (pAM35k_f)

```

                                     CGAATTGATCGGCAAGCCAG   pAM35k_f n = 20 Tm = 56.9°C
gctttgattgccgataactatagctggaatgaatacggcaggtggcaaggctgcctcgaattgatcggaagccagaggcgcgatcgccccctggctgtgatgagacttaga
cgcc
5950      5960      5970      5980      5990      6010      6020      6030      6040      6050
```

pAM35k_f: CGA ATT GAT CGG CAA GCC AG

Reverse primer (pAM35k_r)

```

                                     {CGAATTGATCGGCAAGCCAG}   pAM35k_r n = 20 Tm = 56.9°C
gctttgattgccgataactatagctggaatgaatacggcaggtggcaaggctgcctcgaattgatcggaagccagaggcgcgatcgccccctggctgtgatgagacttaga
cgcc
5950      5960      5970      5980      5990      6010      6020      6030      6040      6050
```

{CGAATTGATCGGCAAGCCAG}

Complement: CTG GCT TGC CGA TCA ATT CG

pAM35k_r: CTG GCT TGC CGA TCA ATT CG

In addition to splitting the pAM backbone in half, it was also split into three parts, with the primers designed as shown below.

Forward primer (pAM27k_f)

```

                                     GCAGGTACATGGGCCAGCAT   pAM27k_f n = 20 Tm = 60.5°C
gccaccagctctcgtagctgctgctgtagctggcaggtacatgggccagcatgggtcgggcgatcgcgacctcatcatcagaagctggagcactgacgatcg
g
5160      5170      5180      5190      5200      5210      5220      5230      5240      5250
```

pAM27k_f: GCA GGT ACA TGG GCC AGC AT

Reverse primer (pAM27k_r)

```

                                     {GCAGGTACATGGGCCAGCAT}   pAM27k_r n = 20 Tm = 60.5°C
gccaccagctctcgtagctgctgctgtagctggcaggtacatgggccagcatgggtcgggcgatcgcgacctcatcatcagaagctggagcactgacgatcg
g
5160      5170      5180      5190      5200      5210      5220      5230      5240      5250
```

pAM27k_r: ATG CTG GCC CAT GTA CCT GC

pAM-23K_f GCTGATATTCCTGCTCT
COMPLEMENT 5'- AGA GCA GCG GGA ATA TCA GC -3'
pAM-23k r AGAGCAGCGGGAATATCAGC

pAM-48k_r AGCTCTTGATCCGGCAAACA
COMPLEMENT 5'- TGT TTG CCG GAT CAA GAG CT -3'
pAM-48k f TGTTTGCCGGATCAAGAGCT

Forward primer (pAM54k_f)

GAGTGTGAGCTGATACCGCTCGCC pAM27k_f n = 22 T_m = 61.9°C

tgattctgtggataaacgctattaccgcctttgagtgtgagctgataccgctcgccgcagccgaacgacccgagcgcagcgaggt < 160

90 100 110 120 130 140 150

pAM54k_f: GAG TGA GCT GAT ACC GCT CGC C

Reverse primer (pAM54k_r)

{GAGTGTGAGCTGATACCGCTCGCC} pAM27k_f n = 22 T_m = 61.9°C

tgattctgtggataaacgctattaccgcctttgagtgtgagctgataccgctcgccgcagccgaacgacccgagcgcagcgaggt < 160

90 100 110 120 130 140 150

pAM54k_r: GGC GAG CGG TAT CAG CTC ACT C

For Gibson primers, fragments of the two different plasmids being joined must be combined into one primer that will overlap regions of the backbone and product vectors. The basic requirements here are that $n \geq 20$ and $T_m \geq 50^\circ\text{C}$. For identification purposes, these primers will be labeled in an A-B_d fashion, where A is the whole fragment from one plasmid, B is the plasmid to which incomplete fragments of A are attached to the other whole fragment, and d is the direction (Figure 3).




Figure 3. Schematic description of Gibson primer structure and naming convention.

PET-pAM primers

Forward primer

```

catctcagaagaggatctgtaa n = 22 Tm = 51.1°C
cgGCCAATTgcagtgaacaaaaactcatctcagaagaggatctgtaa tctagataggtaatctctgcttaaaagcacagaat < 1310
1230      1240      1250      1260      1270      1280      1290      1300

aagggtgggCGCGCGcacc n = 20 Tm = 70.7°C
aatTgtacaaataa aagggtgggCGCGCGcaccagctttcttGtacaaactcggccctcaggagaaggccatcctgacg < 2560
2490      2500      2510      2520      2530      2540      2550

catctcagaagaggatctgtaa aagggtgggCGCGCGcacc
...cgGCCAATTgcagtgaacaaaaactcatctcagaagaggatctgtaa aagggtgggCGCGCGcaccagctttcttGtacaaactcggccctcaggagaag...

```

PET-pAM_f: CAT CTC AGA AGA GGA TCT GTA AAA GGG TGG GCG CGC CGA CCC

Reverse primer

```

{atggcaacaacgacctatt} n = 20 Tm = 52.4°C
catgaatcacaaagtgatggcaacaacgacctattccaagccagtcgCGCGgtttttctggcacagctgggggggctca < 320
250      260      270      280      290      300      310

{ttaagaaggagcccttcacc} n = 20 Tm = 54.3°C
ccTgtttaactttaagaaggagcccttcaccatgtctaaaggTgaagaattattcaactggTgttGtcccaattttggttgaattagatggTgat < 1840
1750      1760      1770      1780      1790      1800      1810      1820      1830

{ttaagaaggagcccttcaccatggcaacaacgacctatt}
...ccTgtttaactttaagaaggagcccttcaccatggcaacaacgacctattccaagccagtcgCGCGgtttttctggcacagctgggggggctca...

```

PET-pAM_r: AAT AGG TCG TTG TTT GCC ATG GTG AAG GGC TCC TTC TTA A

This was not quite done right, as all that was necessary was to take the complement of the PET-pAM_f primer to get the PET-pAM_r primer. In addition, an extra four base pairs were added to the 5' end of the original sequence to prevent primer-dimer potential, and then the complement was taken.

PET-pAM_r: GGG TCG GCG CGC CCA CCC TTT TAC AGA TCC TCT TCT GAG ATG AGT T

pAM-PET primers

Forward primer

atggcaacaacgacctatt n = 20 T_m = 52.4°C
 catgaatcacaaagtatggcaacaacgacctattccaagccagtcgcgcgcttttctggcacagctgggggggctca < 320
 250 260 270 280 290 300 310
 ttaagaaggagcccttcacc n = 20 T_m = 54.3°C
 ccttgtttaactttaagaaggagcccttcaccatgtctaaagggtgaagaattattcaactggtgtgtccaattttggttgaattagatgggtgat < 1840
 1750 1760 1770 1780 1790 1800 1810 1820 1830
 ttaagaaggagcccttcaccatggcaacaacgacctatt
 ...ccttgtttaactttaagaaggagcccttcaccatggcaacaacgacctattccaagccagtcgcgcgcttttctggcacagctgggggggctca...
pAM-PET_f: TTA AGA AGG AGC CCT TCA CCA TGG CAA ACA ACG ACC TAT T

Original reverse primer

{catctcagaagaggatctgtaa} n = 22 T_m = 51.1°C
 cggccaattgcagtgaaacaaaactcatctcagaagaggatctgtaattctagataggtaatctctgctttaaagcacagaat < 1310
 1230 1240 1250 1260 1270 1280 1290 1300
 {aagggtgggcgcgcgaccc} n = 20 T_m = 70.7°C
 aattgtacaaataaagggtgggcgcgcgacccagctttctgtacaaactcggccctcaggagaaggccatcctgacg < 2560
 2490 2500 2510 2520 2530 2540 2550
 {catctcagaagaggatctgtaa} aagggtgggcgcgcgaccc}
 ...cggccaattgcagtgaaacaaaactcatctcagaagaggatctgtaaaagggtgggcgcgcgacccagctttctgtacaaactcggccctcaggagaag...
pAM-PET_r: GGG TCG GCG CGC CCA CCC TTT TAC AGA TCC TCT TCT GAG ATG

Again, the complement of the pAM-PET_f primer should have been taken for the pAM-PET_r primer.
 However, to ensure minimal primer-dimer potential an extra three base pairs were added to the 5' end
 of the original sequence before taking the complement
pAM-PET_r: AAT AGG TCG TTG TTT GCC ATG GTG AAG GGC TCC TTC TTA AAG

Internal PETase primers

Strategy

We do initial screening with the sequencing primers that flank the gene insertion site in the vector backbone. Every plasmid backbone has a _f and _r sequencing primer which is used to sequencing inserted genes, and can also be used to PCR the interest up. When screening by PCR we look for a specific size (size of insert plus flanking sequence) . However this does not provide specific amplification due to the gene, and might be misleading if there is an insert of similar size.

To confirm OUR gene is inserted, we want a specific primer, which will ONLY amplify our gene of interest. This will be a lot more diagnostic. We can do this by making primers specific to the gene, and using one of these as the _f or _r primers in conjunction with the sequencing _r and _f primers.

PETec-N_r1

	Sequence (5'→3')	Template strand	Length	Start	Stop	T _m	GC%	Self complementarity	
Forward primer	CATCGGTGATGTCGGCGATATAGG	Self 3' complementarity	Plus	24	286	309	63.01	54.17	7.00
Reverse primer	CGGACCCCACTTAATGC		Minus	20	777	758	61.03	60.00	4.00
Product length	492								2.00

Primer pair 6

PETec-N_r2

	Sequence (5'→3')	Template strand	Length	Start	Stop	Tm	GC%	Self complementarity		
	Self 3' complementarity									
Forward primer	CATCGGTGATGTCGGCGATATAGG	Plus	24	286	309	63.01	54.17	7.00	2.00	
Reverse primer	GTGTACCCGGGGACGATTG	Minus	19	743	725	60.15	63.16	8.00	2.00	
Product length	458									

PETse-int-C_f

Primer pair 1										
Sequence (5'→3')	Template strand	Length	Start	Stop	Tm	GC%	Self			
complementarity	Self 3' complementarity									
Forward primer	CAGTTCGGTCACAGTTCCCA	Plus	20	761	780	59.89	55.00	3.00	0.00	
Reverse primer	CTTCAAAAAGGCCATCCGTC	Minus	20	1173	1154	57.36	50.00	4.00	1.00	
Product length	413									

PETse-int-N_r

Primer pair 1										
Sequence (5'→3')	Template strand	Length	Start	Stop	Tm	GC%	Self			
complementarity	Self 3' complementarity									
Forward primer	TCATCGGCTCGTATAATGGT	Plus	20	1	20	56.22	45.00	4.00	0.00	
Reverse primer	CACAATTGCGATAGCACCGA	Minus	20	440	421	58.99	50.00	6.00	0.00	
Product length	440									

pETec_C-f3

agtttctggaaattaacggcggg n=23 T_m =57.2°C
 aacagtttctggaaattaacggcggtagccactcttgtgccaaacttgggaacagcaaccaggcactgatcgga
 5880 5890 5900 5910

pETec_C-f4

cttcaaccaacttcagcagtg n=22 T_m =55.7°C
 cgccatgggactcttcaaccaacttcagcagtggtaccgtgccgacgtgatt

GFP-N_r

tgttttctagataccagatcatatgaaacaaca

{tgttttctagataccagatcatatgaaacaaca} n=35 T_m=57.9 C
 ttcaatgtttttctagataccagatcatatgaaacaacatgactttttcaagtctgcatgcc
 1910 1920 1930 1940 1950 1960 1970 1980 1990
 complement: TGT TGT TTC ATA TGA TCT GGG TAT CTA GAA AAA CA

Primer (PCR, Gibson)	Sequence
Fragmentation	
pCS-PET_f	ATG GCA AAC AAC GAC CTA TTC C
pCS-PET_r	TTA CAG ATC CTC TTC TGA GAT GAG TT
pAM_f	AAG GGT GGG CGC GCC GAC
pAM_f2	AAGGGTGGGCGCGCCGACCCAGCTTT

pAM_r	GGT GAA GGG CTC CTT CTT AAA G
PET-pAM_f	CAT CTC AGA AGA GGA TCT GTA AAA GGG TGG GCG CGC CGA CCC
PET-pAM_r	AAT AGG TCG TTG TTT GCC ATG GTG AAG GGC TCC TTC TTA A
pAM-PET_f	TTA AGA AGG AGC CCT TCA CCA TGG CAA ACA ACG ACC TAT T
pAM-PET_r	GGG TCG GCG CGC CCA CCC TTT TAC AGA TCC TCT TCT GAG ATG
pAM35k_f	CGA ATT GAT CGG CAA GCC AG
pAM35k_r	CTG GCT TGC CGA TCA ATT CG
pAM27k_f	GCA GGT ACA TGG GCC AGC AT
pAM27k_r	ATG CTG GCC CAT GTA CCT GC
pAM23k_f	GCTGATATTCCCGCTGCTCT
pAM23k_r	AGAGCAGCGGGAATATCAGC
pAM48k_f	TGTTTGCCGGATCAAGAGCT
pAM48k_r	AGCTCTTGATCCGGCAAACA
pAM54k_f	GAG TGA GCT GAT ACC GCT CGC C
pAM54k_r	GGC GAG CGG TAT CAG CTC ACT C
<i>Screening</i>	
PETec-N_r1	CGGACCCCACTTAATGC
PETec-N_r2	GTGTACCCGGGACGATTG
PETse-int-N_r	CACAATTGCGATAGCACCGA
PETse-int-C_f	CAGTTCCGTCACAGTTCCCA
pETec_C-f3	agtttctggaaattaacggcgggt
pETec_C-f4	cttcaaccaacttcagcagtg
GFP-N_r	TGT TGT TTC ATA TGA TCT GGG TAT CTA GAA AAA CA

(above) Master table of primers used in experiments.

(below) Master table of fragment names with corresponding primers and size. Highlighted are the fragments used to put together pAM-PET.

Name	Template	Forward Primer	Reverse Primer	Size (bp)
PET	pCS_PET	pAM_PET_f	PET_pAM_r	1034
Frag A	pAM	pAM_F2	pAM_35k_r	3521
Frag B	pAM	pAM_35k_f	pAM_r	3554
Frag C	pAM	pAM_F2	pAM_23k_r	2468
Frag D	pAM	pAM_23k_f	pAM_48k_r	2330
Frag D'	pAM	pAM35k_f	pAM48k_r	1278
Frag E	pAM	pAM48k_f	pAM_r	2295

