

Week of June 10, 2019 : Primers were designed and ordered. Cyanobacteria was received.
First attempts at PCR of pAM and PET.

Monday 06/10	1) Designed primers (DK & SM)	Next Steps: Check and order primers						
Tuesday 06/11	1) Received cyanobacteria(NM) 2) ordered primers (NM)	Next Steps: Need to gather protocol for keeping bacteria alive						
Thursday 06/13	1) Received primers (NM)	Next Steps: PCR up fragments (3 fragments)						
Friday 06/14	1) PCR of PET (SH)	Next Steps: 1) PCR frag of pAM 2) run a gel of PCR'd PET						
Saturday 06/15	1) Prepared gel recipe (SM) 2) PCR of pAM (fragment A and B, see details in table 1 below) (SM & SH) 3) ran a gel of PCR'd PET from friday--unsuccessful :(Next Steps: 1) another PCR of PET 2) run a gel of PCR'd Frag A and B						
	Notes: PET PCR was unsuccessful							
	<table><tr><td colspan="2">Table 1</td></tr><tr><td>Frag A</td><td>pAM_r and pAM35k_f</td></tr><tr><td>Frag B</td><td>pAM_f and pAM35k_r</td></tr></table>		Table 1		Frag A	pAM_r and pAM35k_f	Frag B	pAM_f and pAM35k_r
	Table 1							
Frag A	pAM_r and pAM35k_f							
Frag B	pAM_f and pAM35k_r							
Sunday 06/16	1) Attempted to make the gel fragments. We ran 2 tests. 2) This will result in us needing new primers because the current ones are not helping us run proper fragments	Next Steps:						

Week of June 17, 2019 : Cyanobacteria subcultured. PET fragment completed and ready for Gibson assembly. New pCS primers were designed and ordered.

Monday 06/17	AE, DK, SH, SM - 3x PCR of PET	Next Steps: Gel of pAM frag A and B, SM&SK's practice PCR, 3x PET
	Notes: AE did a 1µL Q5 Poly and DK & SH did a .5µL Q5	
Tuesday 06/18	(9:30 - 11:40 AM) SH, AE - Gel electrophoresis of PET(x3), Fragments A and B: PET successful, Fragments A and B unsuccessful (1-1:45 PM) Make a 100 microliter PCR of PET that worked for purification - DK & AE, MD observed	Next Steps 1) Retry Fragments A and B 2) purify PET frag (100 ul new +2*20 ug initial PCR) 3) subculture cyanobacteria with BG11 medium
	Notes: - Most likely need to design smaller A and B fragments to get it to work (3 2200-2500 frags) - SM and EK went to lab to do subculturing, but couldn't find the autoclaved BG11	
Wednesday 06/19	1) Purified PET fragments and checked spectroscopy for success -SM & AE 2) made BG11 plates - EK	Next Steps 1) retry frags A and B 2) design new pCS primers
	Notes: - Our PET fragment is complete and ready for the gibson assembly as soon as we finish our pCS fragments - NM said that when he ran his own personal PCRs of pCS and PET, two worked and one did not - Checked off "purify and validate PET frag" from card "PETase into vector"	
Friday 06/21	PCR for Frag A and B (2 samples each) - SH	Next Steps Gel for Frag A and B 2nd try
	Notes: Samples are labeled A1, A2, B1, B2 along with SH initials	
Sunday 06/23	Gel for Frag A and B 2nd try (unsuccessful) - SH, AE, JH	Next Steps Order new primers, repeat PCR and gel for smaller size

Week of June 24, 2019: Ran PCR of fragments A-E and PET (control).

Monday 06/24	Turbidity reading for cyanobacteria	Next Steps:								
Tuesday 06/25	PCR for PET (control) and Fragments A and B- SH	Next Steps								
	Notes: 2 samples of each, labeled 1 and 2	Gel; if does not work must order new primers -AE								
Wednesday 06/26	1) gel of PET, Frag A and B - AE JH 2) PCR of Frag C, D, and E and PET -SM JH 3) Some cyano preservation - EK	Next Steps: 1) Gel of Frag C, D, E and PET 2) cyano preservation								
	Notes: - Could not find methanol for cyanopres. - EK took turbidity measurements of culture - see table 2 for frag descriptions - 2 samples of each, labeled 1 and 2									
	<table><tr><td colspan="2">Table 2</td></tr><tr><td>Frag C</td><td>pAM_f and 23k_r</td></tr><tr><td>Frag D</td><td>23k_f and 48k_r</td></tr><tr><td>Frag E</td><td>48k_f and pAM_r</td></tr></table>		Table 2		Frag C	pAM_f and 23k_r	Frag D	23k_f and 48k_r	Frag E	48k_f and pAM_r
	Table 2									
Frag C	pAM_f and 23k_r									
Frag D	23k_f and 48k_r									
Frag E	48k_f and pAM_r									
Thursday 06/27	1) Gel of fragments C, D, and E with PET control - AE 2) took cyano light and spect readings -SM	Next Steps: 1) PCR of fragments/purification (?) 2) cryopreservation of cyanobacteria								
	Notes: - see cloning document - gave EK cyano light/spect data									
Saturday 06/29	1) took cyano light and spect readings -SM/SH	Next Steps:								

Week of July 1, 2019 : Ran PCR, gel, and amplification of fragments A, D' and E.

Monday 07/01	PCR of frags E and D' -AE	Next Steps
Tuesday 07/02	Gel of fragments A, D' and E twice -AE	Next Steps : 1) Split cyano culture 2) frag A amplification and purification 3) PCR frag D' and E
	Notes: -Gel unsuccessful both times, retry PCR and gels	
Wednesday 07/03	1) PCR of D' and E + amplification of A -AE + SM	Next Steps : 1) Gel of PCR frag
Thursday 07/04	-gel of frags D', E, A - SM	Next Steps: Amplify frags A, D' and E
Saturday 07/06	1) Amplify Frag A D' and E -SH SM	Next Steps: 1) Run gel of frag A D' and E

Week of July 8, 2019: Successful PCR confirmed by gel for fragments A and D'. PCR unsuccessful for fragment E.

Tuesday 07/09	1) Retest frag A, amplify frag E -AE 2) Gel of frag A and E - DK SM	Next Steps : 1) figure out the proper primers to use for PCR 2) redo PCR
	Notes: - Its highly likely the amplification for E did not work, must be reamplified - frag A ran with program Q5_55_120 - gel of Frag A and E didn't work - If A DOES works, amplify A, D' and E using the new cycle for A. if A DOES NOT work, amplify C, D and E -DONT DO THE ABOVE JUST YET, need to really figure out what primers we need to use for each frag to put it all together	
Wednesday	1) Split liquid cyano culture - SH, EK	Next Steps:

07/10	2) New streak plates done to get individual colonies - SH, EK 3) PCR of B and amplification of A (with high fidelity) -SM SH JH 4) PCR of B and amp of A D' and E (with hot start high fidelity) -SM	1. Make more Bg11 media (need ferric ammonium citrate for this) 2. If individual colonies acquired, can use control plasmid to see if there are problems with uptake into cyanobacteria 3. Gel of PCR set 1 with high fidelity and set 2 with hot start
	Notes: - new streak plates are COLOR CODED according to where they sit on the cabinet - RED should stay closest to the window - BROWN is on the edge of the yellow cabinet - PURPLE is set up on the white block of foam in the middle of the cabinet - See bottom of page for primers used	
Thursday 07/11	1) Gel of PCR of B (with hf and hs) & amplification PCR of A(hf and hs) D' E (hs) -SM 2) PCR of A D' E -SM	Next Steps : 1) Gel of A D' E
	Notes: - Gel unsuccessful, again - Redid A D' E with menhart, if doesn't work we'll need to figure out what's happening	
Friday 07/12	1) Gel of A, D', and E - SH	Next Steps : 1) amplify A D' (if not already) 2) purify frag A D' 3) rerun gel of E
	Notes: - NM is rerunning frag E tom am (hopefully) so we know what's going on with that - successful gel of A and D'	

Week of July 15, 2019: Primers for fragment E redesigned and ordered. Preliminary conjugation protocol developed.

Monday 07/15	Redo gel of fragments A, D', and E samples from 07/11 - SH, AE, SE, DK (E still unsuccessful)	Next Steps : Try a gel previous samples of E that worked Amplification of A and D'
	Notes: If other sample of E works and bp length checks out then we are on our way with fragments	
Tuesday	-Gel of E samples from date - SH -Amplification of A, D' with primers pAM35k_r and pAM_f for	Next Steps : -gel for A and

07/16	frag A and pAM35k_r and pAM48k_f for D' but ran on the wrong cycle. -AE+PC -Gel for frags A, D' and E - AE -amplification of A and D' with the same primers and used the _Q5_57_120 cycle -AE	D' -Redesign primers for frag E
	Notes: All fragments unsuccessful for the first PCR	
Wednesday 07/17	Gel for fragment and of A and D' samples.-PC	Next Steps: Test about the template.
	Notes: Fragments of A and D' still unsuccessful for the PCR.	
Thursday 07/18	-Made a new template DNA and redo the Gels for D',but no result.-PC -Made a new blank gel and test it with standard, only a weak result. -PC	Next Steps -Try a new gel
	Notes: New template with no result, and weak result for standard.	
Sunday 07/21	- started a protocol for conjugation & made new primers for E -SM, SH	Next Steps: - order new primers - confirm protocol
	Notes: - need e. Coli to try conjugation (where's it at?)	

Week of July 22, 2019: First successful PCR for fragments A and B. Amplification was run with taq and Q5 for both fragments.

Monday 07/22	PCR of frag B -SM SH PC	Next Steps: gel of frag B
Tuesday 07/23	1) ran gel of frag b -PC	Next Steps:
	Notes: - Frag B= unsuccessful - rerun with NM	
Wednesday 07/24	1) PCR of multiple Frag -PC NM 2) started 15 ml culture of pet-PET -SM AE PC NM 3) Gel of PCR -SM AE PC	Next Steps: 1) Purify pet-PET culture 2) NM need to review gel to

	<div data-bbox="540 212 787 457" data-label="Image"> </div> <table data-bbox="795 245 1278 457"> <tr> <th>lane</th><th>P_f</th><th>P_r</th><th>expected (approx)</th><th>class name??</th><th>result</th></tr> <tr><td>1</td><td>1kb</td><td>-</td><td></td><td>-</td><td>-</td></tr> <tr><td>2</td><td>100bp</td><td>-</td><td></td><td>-</td><td>-</td></tr> <tr><td>3</td><td>63</td><td>pAM_r</td><td>?</td><td>E?</td><td>ok</td></tr> <tr><td>4</td><td>35</td><td>63</td><td>?</td><td>E?</td><td>no good, two bands</td></tr> <tr><td>5</td><td>F2</td><td>35</td><td>3500</td><td>A</td><td>ok</td></tr> <tr><td>6</td><td>F2</td><td>23</td><td>2300</td><td>C</td><td>maybe ok, but has two bands</td></tr> <tr><td>7?</td><td>35</td><td>pAM_r</td><td>3500</td><td>B</td><td>ok</td></tr> <tr><td>8?</td><td>48</td><td>pAM_r</td><td>2300</td><td>E?</td><td>ok</td></tr> </table> <p data-bbox="540 504 625 529">Notes:</p> <ul data-bbox="587 537 1278 749" style="list-style-type: none"> - FINALLY HAVE SUCCESSFUL FRAGS!!!! (A & B) <ul style="list-style-type: none"> - For some PCR used Q5, for others did taq - Still working on whether to use taq or Q5 - Transforming pET-PET into competent E. Coli to asses PETase 	lane	P_f	P_r	expected (approx)	class name??	result	1	1kb	-		-	-	2	100bp	-		-	-	3	63	pAM_r	?	E?	ok	4	35	63	?	E?	no good, two bands	5	F2	35	3500	A	ok	6	F2	23	2300	C	maybe ok, but has two bands	7?	35	pAM_r	3500	B	ok	8?	48	pAM_r	2300	E?	ok
lane	P_f	P_r	expected (approx)	class name??	result																																																		
1	1kb	-		-	-																																																		
2	100bp	-		-	-																																																		
3	63	pAM_r	?	E?	ok																																																		
4	35	63	?	E?	no good, two bands																																																		
5	F2	35	3500	A	ok																																																		
6	F2	23	2300	C	maybe ok, but has two bands																																																		
7?	35	pAM_r	3500	B	ok																																																		
8?	48	pAM_r	2300	E?	ok																																																		

determine next steps

- Large scale PCR of working frag== FRAG A & B
- Check on gel
- Purify

3) transform pET-PET into competent T7 cell

<p>Thursday 07/25</p>	<p>1) Amplification for Frag A & B with different enzyme both(Q5 and taq) -AE PC</p> <table data-bbox="540 896 1278 1180"> <tr> <th>Fragment :</th><th>w/taq</th><th>w/Q5</th></tr> <tr> <td>A</td><td>pAM_F2 + 35k_r</td><td>pAM_F2 + 35k_r</td></tr> <tr> <td>B</td><td>35k_f + pAM_r</td><td>35k_f + pAM_r</td></tr> </table> <p>2) Purified pet-PET culture -AE PC</p> <p data-bbox="540 1285 625 1312">Notes:</p> <ul data-bbox="540 1358 1221 1428" style="list-style-type: none"> - Documented the "PCR Purification" on google drive - did amplification of BOTH frags with taq AND Q5 	Fragment :	w/taq	w/Q5	A	pAM_F2 + 35k_r	pAM_F2 + 35k_r	B	35k_f + pAM_r	35k_f + pAM_r
Fragment :	w/taq	w/Q5								
A	pAM_F2 + 35k_r	pAM_F2 + 35k_r								
B	35k_f + pAM_r	35k_f + pAM_r								

Next Steps :

- Gels for the frags
- Replicas plate NEB3013


<p>Friday 07/26</p>	<p>1) Gels of Frag A , B, PET(taq)</p> <table data-bbox="540 1575 1196 1785"> <tr> <th></th><th>lane</th><th>P_f</th><th>P_r</th><th>expected (approx)</th><th>class name??</th><th>result</th></tr> <tr><td></td><td>1</td><td>-</td><td>-</td><td></td><td>-</td><td>-</td></tr> <tr><td></td><td>2</td><td>100bp</td><td>-</td><td></td><td>-</td><td>-</td></tr> <tr><td></td><td>3</td><td>pAM_f</td><td>pAM_r</td><td>1034</td><td>PET</td><td>ok</td></tr> <tr><td>taq</td><td>4</td><td>F2</td><td>35</td><td>3521</td><td>A</td><td>maybe ok, the size might be too big</td></tr> <tr><td>taq</td><td>5</td><td>35</td><td>pAM_R</td><td>3554</td><td>B</td><td>no</td></tr> <tr><td>Q5</td><td>6</td><td>F2</td><td>35</td><td>3521</td><td>A</td><td>maybe ok, but has two bands</td></tr> <tr><td>Q5</td><td>7?</td><td>35</td><td>pAM_r</td><td>3554</td><td>B</td><td>ok</td></tr> </table> <p data-bbox="540 1827 625 1854">Notes:</p>		lane	P_f	P_r	expected (approx)	class name??	result		1	-	-		-	-		2	100bp	-		-	-		3	pAM_f	pAM_r	1034	PET	ok	taq	4	F2	35	3521	A	maybe ok, the size might be too big	taq	5	35	pAM_R	3554	B	no	Q5	6	F2	35	3521	A	maybe ok, but has two bands	Q5	7?	35	pAM_r	3554	B	ok
	lane	P_f	P_r	expected (approx)	class name??	result																																																			
	1	-	-		-	-																																																			
	2	100bp	-		-	-																																																			
	3	pAM_f	pAM_r	1034	PET	ok																																																			
taq	4	F2	35	3521	A	maybe ok, the size might be too big																																																			
taq	5	35	pAM_R	3554	B	no																																																			
Q5	6	F2	35	3521	A	maybe ok, but has two bands																																																			
Q5	7?	35	pAM_r	3554	B	ok																																																			

Next Steps:

- Replicas plate NEB3013
- Purify the B(taq) and both A,B(Q5) - NM DECISION
- Rerun the gel of A(pAM_f2 + 35k_r, taq) (1 ul) with 1kb and 1ul of purified samples. - NM

	<p>-Most of the gels are good, but A(pAM_f2 + 35k_r, Q5) shows 2 bands B(taq) no band A(pAM_f2 + 35k_r, taq) may be too big -No colonies shows off on the plate NEB3013, need wait more 24 hours and see.</p>	
--	---	--

Week of July 29, 2019: Purification of fragments QA, QB and taq B. Gibson reaction of the primer A, B , pet-PET documented. Started liquid culture of pet_PET/3013.

<p>Tuesday 07/30</p>	<p>1) Reran the gel of A(pAM_f2 + 35k_r, taq) (1 ul) with 1kb and 1ul of purified samples -SM AE PC</p> <table><thead><tr><th></th><th>lane</th><th>P_f</th><th>P_r</th><th>expected (approx)</th><th>class name??</th><th>result</th><th></th></tr></thead><tbody><tr><td></td><td>1</td><td>1kb</td><td>-</td><td></td><td></td><td>-</td><td></td></tr><tr><td></td><td>2</td><td>100bp</td><td>-</td><td></td><td></td><td>-</td><td></td></tr><tr><td>taqA (unpu</td><td>3</td><td>F2</td><td>35</td><td>3500</td><td>A</td><td></td><td></td></tr><tr><td>taqB</td><td>4</td><td>35</td><td>pAM_R</td><td>3550</td><td>B</td><td></td><td></td></tr><tr><td>QA</td><td>5</td><td>F2</td><td>35</td><td>3500</td><td>A</td><td></td><td></td></tr><tr><td></td><td>6</td><td>blank</td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>QB</td><td>7</td><td>35</td><td>pAM_r</td><td>3550</td><td>B</td><td></td><td></td></tr></tbody></table> <p>NOTES: -the gel was ripped in lanes 7 and 6, which could explain why the band is so light - there is no band for frag B w/ taq (lane 4) and frag A w/ Q5 has double bands - all bands seem to be similar in size which corresponds to their expected sizes</p> 		lane	P_f	P_r	expected (approx)	class name??	result			1	1kb	-			-			2	100bp	-			-		taqA (unpu	3	F2	35	3500	A			taqB	4	35	pAM_R	3550	B			QA	5	F2	35	3500	A				6	blank						QB	7	35	pAM_r	3550	B			<p>To Do:</p> <ul style="list-style-type: none">- need info/assistance with NEB3013- need info/assistance on concentrations for GFP conjugation media- go over gel from 0730
	lane	P_f	P_r	expected (approx)	class name??	result																																																												
	1	1kb	-			-																																																												
	2	100bp	-			-																																																												
taqA (unpu	3	F2	35	3500	A																																																													
taqB	4	35	pAM_R	3550	B																																																													
QA	5	F2	35	3500	A																																																													
	6	blank																																																																
QB	7	35	pAM_r	3550	B																																																													
<p>Thursday 08/01</p>	<p>1) spec of purified frags (QB QA and taq B) 2) pcr of taqB -purify taq A and Q5B -PC</p> <p>Notes: - using taq A and Q5B with pet-PET for plasmid</p>	<p>Next Steps :</p> <ol style="list-style-type: none">1) gel of purified fragments2) redo competency for pet_PETase3) try conjugation with e. Coli4) rewrite “intro to project” for JH<ul style="list-style-type: none">- Powerpoint of stuff5) gather docs/pics for website																																																																

Friday 08/02	-gel of purified frags(Q5B taqA) -PC -Gibson reaction of the primer A,B,pet-PET -PC	Next Steps :
	-Documented Gibson reaction setup in the lab manual.	
Saturday 08/03	1) replicata plated pet_PET/3013 -SM	Next Steps : 1) need to grow pet_PET/3013 culture in liquid media 2) may need to make 100/200 ul of Frag
	Notes: - strong colonies chosen for replication; placed replica plate in incubator - gibson plate put back into incubator; no colonies at 1100a, had strong colonies on the control plate	
Sunday 08/04	1) started liquid culture of pet_PET/3013 -SM	Next Steps : 1) need 200ul rxn of frag A and B
	Notes: - strong colonies of pet_PET/3013 on LBa plate - liq culture of pet_PET/3013 in shaker in 311 - LBa plate w/ pet_PET/3013 in fridge with LBc plate - one small, questionable colony on gibson plate (pAM_PET/NEB5a on LBc) - pAM_PET/NEB5a in incubator, control on lab bench	

Week of August 5, 2019: PCR and purification of fragments A and B with Q5. Gibson reaction ran and transformed.

Monday 08/05	1) prepped gel -AE 2) PCR: 150ul taqA, 100ul Q5B 3) plated pAM4788	Next Steps : 1) 20 ul of Q5A (gave funny two bands, refer to pAM primer test doc) with mastermix 2) make mini preps of pet_PET/3013
	Notes: - didn't have enough enough template for 20 ul rxn, made some new plates today (0805), mini prep tomorrow (08	
Tuesday	1) Ran gel for amplified fragments A and B -AE 2) Purified fragments A and B and took a spec	Next Steps :

08/06	reading both -AE 3) Prepared new template for PCR - AE and NM 4) Set up and incubated Gibson reaction -AE SM	1) Mini prep plasmids from over the weekend
	Notes: - start trying to do conjugation between e. Coli and cyano -- make media -stored the remaining gibson in the -20 freezer	
Wednesday 08/07	Notes: - no growth on plates AE SM did with gibson	
Thursday 08/08	1) replated gibson rxn -NM 2) PCR of frag A and B with Q5 mastermix	Next Steps: - mini prep plasmids from weekend (pet_PET/3013)
	Notes: - SM AE didn't transfer bacteria onto plate correctly, NM took the stored gibson bacteria and replated (who knows if the bacteria will grow after freezing)	
Friday 08/09	1) PCR of frag A and B with Q5 -NM 2) Purification of PCR -SM 3) Gibson -SM	Next Steps : 1) transformation from gibson
	Notes: - Plate with gibson from tues 06 didn't take	
Saturday 08/10	1) transformation from gibson -SM	

Week of August 12, 2019: PCR of Gibson reaction (also cultured and plated colonies). No successful Gibson reaction with the correct plasmid.

Monday 08/12	1) PCR of Gibson reaction -AE 2) Cultured colonies from Gibson plate -AE 3) Plated the colonies from the Gibson -AE 4) 2 gels for the PCR -AE	Next Steps : - mini prep plasmids (pet_PET/3013) - website - redo gibson ???
	Notes: - no successful gibson with CORRECT plasmid	

Week of August 19, 2019: Did mini preps of pet_PET and Gibson transformation.

Saturday 08/24	- mini preps of pet_PET(red minicent tubes) and gibson transformation (blue minicent. tubes) -SM JT WJB	Next Steps : - screen pET-PET
-------------------	---	----------------------------------


		and pAM_PET - redo gibson ???
	Notes: - refer back to gel doc located ***** for replica plate colonies for gibson	

Week of August 26, 2019: Re-ran gel of pET_PET to ensure colony growth. Made culture of pET_PET colonies grown from replica plate.

Monday 08/26	Made cyano medium, checked progress of subculture - EK MD	Next Steps :
Tuesday 08/27	- screened mini preps from aug 24 -SM AE	Next Steps : - run a gel of PCR - finalize a protocol for GFP transformation between e.coli and cyano
	Notes - ran: TAQ_57_120 to screen	
Wednesday 08/28	-Assay team touched base with Matt from summer	Next Steps :
Friday 08/30	- re-ran gel of pET_PET to make sure there was a colony -SM AE JT ML	Next Steps : - grow culture of pET_PET colonies from replica plate
	Notes: - may have some colonies with pET_PET, confusion on which so growing up two	
Saturday 08/31	- plated colonies 6,7, and control for pET_PET -SM	Next Steps : - Move pET_PT plates to fridge
	Notes: - plates in incubator, need to be moved to fridge Sun.	
Sunday 09/01	- pET_PET plates in fridge -SM	Next Steps :

Week of September 2, 2019: Made a liquid culture of pET_PET and did mini preps of four cultures (-6, -7, C, and C). Another Gibson reaction was set up and transformation was run. By the end of the week, no transformed colonies were observed.

	Assay Team: - liq culture of pET_PET -SM SS WB RB	Next Steps :
--	---	--------------

Tuesday 09/03		- minipreps of 4 cultures																																	
Wednesday 09/04	Assay : Miniprep of four cultures -SS,CJ,LG	Next Steps : Cloning: Gibson Assay: PCR screen, Protein mini prep																																	
	Notes: minipreps left in incubator overnight,suggested by NM, labeled -6,-7,C,C.																																		
Thursday 09/05	Cloning: - set up gibson rxn -SM JT RB - transformation -SM AE - made LBsm plates (~24) Assay: - made liq cultures of -6, -7, -NOs -SM WB AE	Next Steps : Cloning: Assay: - mini preps of liq cultures (protein screen) - PCR screen																																	
	Notes: Gibson: - "NEBulider DNA assembly master mix" - "NEBuilder positive control" - 2 control tubes (labeled: +)(one SM one RB) - 1 pAM_PET tube (labeled: G) Assay: - unable to do PCR because water wasn't clean enough																																		
Friday 09/06	Cloning: - no observed transformed colonies on plate @ 1410pm-AE Assay: - screen PCR -SM WB - screen PCR gel of EN and EC - ML - Protein mini preps -JC <table><tr><td>lane</td><td>P_f</td><td>P_r</td></tr><tr><td>1</td><td>pet_f</td><td>petec_N_r1</td></tr><tr><td>2</td><td>"</td><td>"</td></tr><tr><td>3</td><td>"</td><td>"</td></tr><tr><td>4</td><td>"</td><td>"</td></tr><tr><td>5</td><td>petec_C_f1</td><td>pet_r</td></tr><tr><td>6</td><td>"</td><td>"</td></tr><tr><td>7</td><td>"</td><td>"</td></tr><tr><td>8</td><td>"</td><td>"</td></tr><tr><td>9</td><td></td><td></td></tr><tr><td>10</td><td></td><td></td></tr></table> 	lane	P_f	P_r	1	pet_f	petec_N_r1	2	"	"	3	"	"	4	"	"	5	petec_C_f1	pet_r	6	"	"	7	"	"	8	"	"	9			10			Next Steps : Cloning: - check for growth later in the day - if growth take out and replica plate. Replica plate goes in incubator and old plate goes in fridge - if no growth, figure out what is wrong with gibson (consult NM)
lane	P_f	P_r																																	
1	pet_f	petec_N_r1																																	
2	"	"																																	
3	"	"																																	
4	"	"																																	
5	petec_C_f1	pet_r																																	
6	"	"																																	
7	"	"																																	
8	"	"																																	
9																																			
10																																			
	Notes: Assay: - questionable mini prep of -7 and -C2																																		



	<ul style="list-style-type: none"> - PCR tubes under "EE" (with pET_f/r) had drastic volume differences didn't run a gel - used primer pET_f rather than pET_f2 	
--	---	--


Week of September 9, 2019: Redid Gibson reaction with LBamp plates for control and performed transformation procedure. No growth was observed on Gibson plates once again. Made cultures of pCS_PET and pAM4788. Ran PCR for PET, A, and B. Also ran PCR for EE screening of -6, -7, C1, C2 from pET_PET and did protein purification.

Monday 09/09	Cloning: - took out neg gibbon plates -SM	Next Steps: Cloning: - redo gibbon one more time
	Notes: Cloning: <ul style="list-style-type: none"> - Had "positive" colonies last time we ran gibbon, why no colonies this time - Used LBsm plates for control, needs to be LBamp 	
Tuesday 09/10	Cloning: <ul style="list-style-type: none"> - redid gibbon. one experimental tube and one positive control -WB MD - Transforming gibbon -SM ML - Made gel for frag a b pet check -SM ML 	Next Steps: <ul style="list-style-type: none"> - Ran a gel of fragments A, B and PET - Check plates in incubator
	Notes: <ul style="list-style-type: none"> - NEED TO ORDER MORE C2987 CELLS FROM NEB - PCR used flat lids instead of rounded, so tube became malformed. Took a lot of work to get product ou 	
Wednesday 09/11	Cloning: <ul style="list-style-type: none"> - no growth on gibbon plates -SM - conformation gel of A B and PET (see gel below) -SM 	Next Steps: Cloning: <ul style="list-style-type: none"> - make frag A and B and PET via PCR - talk to NM tomorrow about what to do with no gibbon growth
	Notes: Cloning: <ul style="list-style-type: none"> - no band for PET, may be the cause of failed gibbon reactions - control gibbon grew, no rxn for pAM_PET 	
Thursday	Cloning:	Next Steps:

09/12	- replated pCS_PET and pAM4788 -SM	Cloning: - PCR of frag A B and PET - if no growth, figure out what is wrong with gibson (consult NM)
	Notes: Cloning: - might need to make new template b/c it's getting old, hence why replated	
Friday 09/13	Cloning: - Ran PCR for PET, A, and B -SS, AE, MD - started pAM4788 and pSC_PET cultures -SM CJ Assay: - PCR for EE screening of -6, -7, C1, C2 from pET_PET -SM - started liq cultures of -6, -7, NOS control, cindel control (for protein prep) -JC	Next Steps: Cloning: - gel of frag A B PET - do PCR purification with carts Assay: - gel of EE screening
	Notes: Cloning: Q5MM with _Q5_57_120 Assay: See below for "EE" PCR parts	
Saturday 09/14	Assay: Ran Protein Purification from Friday's overnight cultures- LG JCT	Next Steps: Run the gel to determine expression
Sunday 09/15	Assay: Vlad set up + ran gel to verify results	


Week of September 16, 2019: Purified fragments A, B and PET. Received a very low A260 value for PET frag, but OK A260 values for fragments A and B. Started new PCR of PET fragment to redo. Transformation of pAM4788 and Gibson into E. coli. Did high quality mini preps of pET_PET for assay.

Monday 09/16	<p>Cloning: - ran gel of Frag A B and PET -SM</p> <p>Assay: Read gel from 15 sept; empty</p>  <table><thead><tr><th></th><th>lane</th><th>P_f</th><th>P_r</th><th>expected</th><th>class</th><th>result</th></tr></thead><tbody><tr><td></td><td>1</td><td>1kb</td><td>-</td><td></td><td></td><td>good</td></tr><tr><td>A</td><td>2</td><td>f2</td><td>35r</td><td>3500</td><td>A</td><td>good</td></tr><tr><td>B</td><td>3</td><td>35f</td><td>pAM_R</td><td>3500</td><td>B</td><td>good</td></tr><tr><td>PET</td><td>4</td><td>pam_petf</td><td>pam_petr</td><td>1034</td><td>PET</td><td>very light, not so great</td></tr><tr><td>A</td><td>5</td><td>f2</td><td>35r</td><td>3500</td><td>A</td><td>not so great, appears to have 3 bands</td></tr><tr><td>B</td><td>6</td><td>35f</td><td>pamr</td><td>3500</td><td>B</td><td>good</td></tr><tr><td>PET</td><td>7</td><td>pam_petf</td><td>pam_petr</td><td>1034</td><td>PET</td><td>light band, but present</td></tr></tbody></table>		lane	P_f	P_r	expected	class	result		1	1kb	-			good	A	2	f2	35r	3500	A	good	B	3	35f	pAM_R	3500	B	good	PET	4	pam_petf	pam_petr	1034	PET	very light, not so great	A	5	f2	35r	3500	A	not so great, appears to have 3 bands	B	6	35f	pamr	3500	B	good	PET	7	pam_petf	pam_petr	1034	PET	light band, but present	<p>Next Steps: Cloning: - PCR purification - gibson/ transformation</p> <p>Assay: - gel of PET “EE” screen</p>
	lane	P_f	P_r	expected	class	result																																																				
	1	1kb	-			good																																																				
A	2	f2	35r	3500	A	good																																																				
B	3	35f	pAM_R	3500	B	good																																																				
PET	4	pam_petf	pam_petr	1034	PET	very light, not so great																																																				
A	5	f2	35r	3500	A	not so great, appears to have 3 bands																																																				
B	6	35f	pamr	3500	B	good																																																				
PET	7	pam_petf	pam_petr	1034	PET	light band, but present																																																				
	<p>Notes: Cloning: - positive bands for A B PET</p>																																																									
Thursday 09/19	<p>Cloning: - purified frag A B PET -SM</p> <p>Assay: -- gel of pET_PET screen of “EE” (pet_f and pet_r) -SM</p> <table><thead><tr><th>Lane</th><th>P_f</th><th>P_r</th><th>size exp</th><th>results</th></tr></thead><tbody><tr><td>1</td><td>pet_f</td><td>pet_r</td><td>1034</td><td>poor</td></tr><tr><td>2</td><td>"</td><td>"</td><td></td><td>poor</td></tr><tr><td>3</td><td>"</td><td>"</td><td></td><td>well defined band</td></tr><tr><td>4</td><td>"</td><td>"</td><td></td><td>well defined band</td></tr><tr><td>5</td><td></td><td></td><td></td><td>-</td></tr><tr><td>6</td><td></td><td></td><td></td><td>-</td></tr></tbody></table> 	Lane	P_f	P_r	size exp	results	1	pet_f	pet_r	1034	poor	2	"	"		poor	3	"	"		well defined band	4	"	"		well defined band	5				-	6				-	<p>Next Steps: Cloning: -- autoclave tips !!!! - purify and make new template of pCS_PET for PET frag - redo PCR of frag PET</p>																					
Lane	P_f	P_r	size exp	results																																																						
1	pet_f	pet_r	1034	poor																																																						
2	"	"		poor																																																						
3	"	"		well defined band																																																						
4	"	"		well defined band																																																						
5				-																																																						
6				-																																																						
	<p>Notes: Cloning: - SUPER low a260 for PET - not great a260 for A and B, but workable</p>																																																									
Friday 09/20	<p>Cloning: - PCR of PET frag -SS WB - Gel of PET frag -SM - started liq culture of pAM4788 (for control transformation into S.e) -SM - Gibson -SM</p>	<p>Next Steps:</p> <p>Cloning: - HQ minipreps of pAM4788 → transformation</p>																																																								

	<div></div> <table><thead><tr><th></th><th>lane</th><th>P_f</th><th>P_r</th><th>expected (approx)</th><th>class name??</th><th>result</th></tr></thead><tbody><tr><td>100</td><td>1</td><td>-</td><td>-</td><td></td><td></td><td>good</td></tr><tr><td>PET frag</td><td>2</td><td>pam_px</td><td>pam_PETf</td><td>1034</td><td></td><td>good!!</td></tr></tbody></table> <p>Assay:</p> <ul style="list-style-type: none">- plated pET_PET in DH5a (og) -SM- liq cultures of pET_PET (for PET screening) -SM <p>Notes:</p> <ul style="list-style-type: none">- need to autoclave pipette tips and make more LBa liq broth		lane	P_f	P_r	expected (approx)	class name??	result	100	1	-	-			good	PET frag	2	pam_px	pam_PETf	1034		good!!	<p>into TURBO</p> <p>→ test conjugation into S.e</p> <p>Assay:</p> <ul style="list-style-type: none">- HQ mini preps of pET_PET→ PCR screen once we have new ec_C_
	lane	P_f	P_r	expected (approx)	class name??	result																	
100	1	-	-			good																	
PET frag	2	pam_px	pam_PETf	1034		good!!																	
Saturday 09/21	<p>Cloning:</p> <ul style="list-style-type: none">- HQ minipreps of pAM4788- transformation of pAM4788 into turbo- transformation of gibson into turbo <p>Assay:</p> <ul style="list-style-type: none">- HQ minipreps of pET_PET	<p>Next Steps:</p> <p>Cloning:</p> <ul style="list-style-type: none">- check transformation plates→ if pos: replica plate and screen pAM4788 and pAM_PET→ if neg: redo transformation of pAM4788 and redo gibson of pAM_PET <p>Assay:</p> <ul style="list-style-type: none">- PCR screen once new C primers come in																					

Week of September 23, 2019: Made mini preps of four successful Gibson colonies and ran PCR screen.

Monday 09/23	<p>Cloning:</p> <ul style="list-style-type: none"> - four isolated colonies on gibson plate (control growth=good) 	<p>Next Steps :</p> <p>Cloning:</p> <ul style="list-style-type: none"> - make minipreps
--------------	--	--

	- growth of pAM4788 on plate	of gibson 4 colonies → run a screening PCR - screen turbo cells for pAM4788 Assay: - order pET C primers - PCR screen																																																
Tuesday 09/24	- mini preps of gibson 4 colonies -SM AE - PCR screen -ML	Next Steps : Cloning: - run gel of PCR screen - order GFP primers to screen pAM4788 - made just need to order Assay: - order pET C primers - made just need to order - PCR screen																																																
	Notes: - no more taq MM-- order more																																																	
Wednesday 09/25	- PCR screen gel -SM	Next Steps: - redo PCR screen																																																
Thursday 09/26	- PCR screen of -1,2,3,4,+ pAM-PET (+pAM4788) -SM <table><thead><tr><th>lane</th><th>P_f</th><th>P_r</th><th>expected (approx)</th><th>class name??</th><th>result</th></tr></thead><tbody><tr><td>1</td><td>pAMse_f</td><td></td><td></td><td></td><td></td></tr><tr><td>2</td><td>"</td><td>pAMse_r</td><td></td><td></td><td></td></tr><tr><td>3</td><td>"</td><td>"</td><td>1200</td><td></td><td>no</td></tr><tr><td>4</td><td>"</td><td>"</td><td>1200</td><td></td><td>no</td></tr><tr><td>5</td><td>"</td><td>"</td><td>1200</td><td></td><td>no</td></tr><tr><td>6</td><td>"</td><td>"</td><td>1200</td><td></td><td>no</td></tr><tr><td>7</td><td>"</td><td>"</td><td>900</td><td></td><td>YES</td></tr></tbody></table> 	lane	P_f	P_r	expected (approx)	class name??	result	1	pAMse_f					2	"	pAMse_r				3	"	"	1200		no	4	"	"	1200		no	5	"	"	1200		no	6	"	"	1200		no	7	"	"	900		YES	
lane	P_f	P_r	expected (approx)	class name??	result																																													
1	pAMse_f																																																	
2	"	pAMse_r																																																
3	"	"	1200		no																																													
4	"	"	1200		no																																													
5	"	"	1200		no																																													
6	"	"	1200		no																																													
7	"	"	900		YES																																													

Week of October 7,2019: Redid Gibson and PCR of fragments A, B, and PET. Did Gibson extraction PCR using magic peg and had 11 colonies grow on the plate.

Wednesday	Redid gibson	Next Steps :
-----------	--------------	--------------

10/09	PCR of frag A <hr/> Notes: <ul style="list-style-type: none"> - Need 35k_f and pAM_PET_r - Frag A= pAM_f2 and 35k_r 	PCR of frag B and PET
Thursday 10/10	- PCR of B and PET -Transformation -Made gel <hr/> Notes Frag B= 35k_f and pAM_r Frag PET= pAM_PET_f and PET_pAM_r	Next Steps : - Run a gel of A B PET -Purify frag
Friday 10/11	- Gibson extraction PCR (magic peg) -Started 3ml liq cultures/ replica plate of 11 gibson colonies -Started 15ml liq culture/plate of #9 colony <hr/> Notes: - Had 11 colonies grow on gibson plate !	Next Steps : - Run a gel of extraction (magic peg) PCR and of A B PET frag
Saturday 10/12	- Ran gel of magic peg PCR -HQ purification of #9 -PCR screen of #9	Next Steps : - Run a gel of A B PET frag Run a gel of PCR screen

Week of October 14, 2019: TFA assay was developed to prove that the degradation of PET nanoparticles tagged with fluorescent molecules could be tracked through fluorimetry.

Monday 10/14	Ran a gel of PCR screen. Ran a gel of Re-PCR screen. Gel of restriction enzyme screen. Restriction enzyme screen.	Next Steps : -Re-PCR screen one with Q5 another with taq
Wednesday 10/16	Developed TFA assay with Dr Menhart.	Next Steps : - Calculate moles of TFA, find appropriate buffer
Thursday 10/17	Calculations done, chose tris. Set up trials with 0/20/45/90% TFA, left in incubator	Next Steps : - Add buffer & check in fluorimeter

